Dry Coat

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English

Product description

Dry Coat is a light curing one-bottle varnish, based on multifunctional methacrylates for coating of surfaces.

Indications

- Coating of temporaries, customized impression trays and ortho dontic plates.
- As a primer/bonding agent for indirect restorations of PEEK, composite, and PMMA

Contra-indications

Not to be used if allergies are known to exist against constituents (especially methylmethacrylate).

Basic safety instructions

- ▶ For dental use only!
- ▶ Keep out of reach of children!

Interactions with other materials

- ▶ Products containing eugenol can affect the curing of materials. Avoid contact with the materials containing eugenol.

Instructions for use

- ▶ Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm².
- ▶ Close the bottle immediately after use.

Recommended use

Surface coating

1. Apply a thin layer of the material without inclusion of air bubbles.
2. Let the material act for 20 s and the cure each element for 20 s with a halogen or LED lamp. Alternatively cure the entire restoration for at least 90 s in the light furnace.

Note: To ensure that no oxygen inhibition layer is formed the material should be light cured within 2 minutes of use.

Bonding agent for restorations

1. Prepare the restoration according to manufacturer's specifications.
2. Apply a thin layer of the material without entrapped air bubbles to the required surfaces. Expose for 20 s with a halogen lamp or ≥ 10 s with an LED lamp (light intensity >1000 mW/cm²).
3. Lute restoration with a suitable adhesive dual- or self-curing composite luting cement according to the manufacturer's instructions.

Storage and shelf life

- ▶ Store in a dry place at 2 to 25 °C/36 to 77 °F!
- ▶ Do not use after the expiration date.

Composition

Multifunctional acrylate, methylmethacrylate, catalyst, stabilizer, additive

Packaging

1 bottle @ 30 ml liquid	REF 829197
-------------------------	------------

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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