Ecosite Bulk Fill

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English

Product description

Ecosite Bulk Fill is a light-curing, radiopaque restorative material developed especially for use in the posterior region.

Indications

- Restorations in the posterior region (Classes I and II, including replacement of individual cusps)
- Restorations of class V (cervical caries, root erosions, wedge-shaped defects)
- Core build-up
- Extended fissure sealing on molars and premolars
- Restorations in deciduous teeth

Contraindications

- Application of the material is contraindicated if dry isolation or the recommended application technique is not possible.
- Do not use directly on exposed pulp or in cases of allergies to any of the ingredients.

Interactions

- Avoid materials containing eugenol, moisture and oily air, as these may inhibit polymerization at the contact area.
- The use of cationic oral rinses, plaque indicators and chlorhexidine may lead to discoloration.

Side effects

There are no known side effects to date. Hypersensitivity to components of the material cannot be ruled out in individual cases.

Basic safety information

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.

Notes for use

- ▶ For optimum polymerization, position the light-curing unit as close to the restorative material as possible and observe the recommendations relating to the light-curing time and the light intensity (see table). Observe the manufacturer's instructions for the light-curing unit used.
- ▶ Follow the manufacturer's instructions for other materials that are used with Ecosite Bulk Fill.
- ▶ If Ecosite Bulk Fill is applied into the patient's mouth straight from the Safetip, the Safetip must only be used on one patient for hygiene reasons.

Recommended use

1. Clean the tooth. Remove all cleaning residue with water spray. It is advisable to use a rubber dam.
2. Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Ensure that the areas close to the pulp are protected.
3. Condition the dental enamel and apply the adhesion primer according to the manufacturer's instructions for the material used. DMG advises using the universal adhesion primer Ecosite Bond.
4. If necessary, apply an underfilling of low-viscosity composite or compomer. We recommend using a matrix.
5. Apply Ecosite Bulk Fill in the cavity with a maximum layer thickness of 5 mm. Apply taller restorations one layer at a time. Light-cure each layer separately.

| | | |
|-----------------|----------------------------|--------------------------|
| Light intensity | > 1,000 mW/cm ² | > 600 mW/cm ² |
| Curing time | 20 s | 40 s |

Note: Additional buccal or lingual light-curing is advisable. Remove metal matrices beforehand.

6. Removal of excess material, contouring and fine processing are done using diamond finishing tools and flexible polishing discs of decreasing grain size, for example. Ecosite Bulk Fill can be given a high shine using polishing paste together with polishing cups or discs.

Disinfection of the Safetip-Dispenser

The Safetip-Dispenser is reusable after disinfection, although it should be replaced if visible signs of wear appear.

Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute). The Safetip-Dispenser is also autoclavable.

Composition

Barium glass in a Bis-GMA-based matrix of dental resins.

Pigments, additives and catalyst. Filler content: 82% by weight = 65% by volume

The range of variation of the inorganic filler particles is between 0.02 and 0.7 µm

Storage and shelf life

- Store at room temperature (2 to 25°C / 36 to 77°F).
- Close syringes / composite tip immediately after use so that the material does not cure in the ambient light.
- Do not use after the expiration date.

Packaging

| | |
|--|------------|
| Starter Kit: | |
| 1 syringe containing 4 g of Light 2 syringes each containing 4 g of Universal 1 bottle @ 5 ml Ecosite Bond | REF 220750 |
| 16 Safetips each containing 0,25 g of Light 32 Safetips each containing 0,25 g of Universal 1 bottle @ 5 ml Ecosite Bond | REF 220751 |

| | |
|--|------------|
| Packaging: 1 syringe containing 4 g paste | |
| Light | REF 220752 |
| Universal | REF 220753 |
| Contrast | REF 220754 |
| Packaging: 16 Safetips each containing 0.25 g paste | |
| Light | REF 220755 |
| Universal | REF 220756 |
| Contrast | REF 220757 |

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