Ionosit-Baseliner

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Product description

Ionosit-Baseliner is a light-curing, radiopaque, compomer base liner. The controlled expansion of Ionosit-Baseliner counteracts the shrinkage of a composite filling.

Indications

- Base liner for composite restorations
- Blocking out of undercuts

Basic safety instructions

- ▶ For dental use only!
- ▶ Keep out of reach of children!

Side effects

There have been no systemic side effects to date. Isolated cases of hypersensitivity or contact allergy to components of the material cannot be excluded.

Instructions for use

- ▶ Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm². Place the light as close as possible to the material.
- ▶ Ionosit-Baseliner adheres to metal instruments. Therefore, clean instruments with a paper towel before the material sets.
- ▶ Ionosit-Baseliner is light-sensitive. Protect the material from premature activation by ambient light and close the syringe immediately after use.
- ▶ Do not remove the oxygen inhibited layer because it is required for bonding to the next layer.
- ▶ The use of a matrix is recommended. Follow the manufacturers' instructions for other materials used.

Recommended use

1. Before preparation, clean the tooth and remove all cleaning residue with water spray. Use of a rubber dam is not necessary if moisture can be controlled in another way. Avoid contamination with saliva during the application.
2. Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Protect pulp in deep cavities with suitable materials.
3. Etch the cavity in accordance with the preferred technique (e.g. Total-Etch). Observe the manufacturer's instructions for the material used.
4. Apply bonding agent (e.g. TECO) in accordance with the manufacturer's instructions.
5. Ionosit-Baseliner can be used both as a conventional base liner (coating the floor of the cavity) and to coat the entire cavity (stress-breaking function). Apply Ionosit-Baseliner with a brush or ball-pointed instrument in a thickness of about 1 mm.
6. Cure for at least 20 seconds with a suitable light unit.
7. Place the composite filling material according to the manufacturer's instructions. DMG recommends the use of EcuSphere.

Storage and shelf life

- ▶ Store in a dry place at room temperature (15 to 25 °C/59 to 77 °F)!
- ▶ Do not use after the expiration date.

Composition

Ionomer glass in a matrix of polymerizable oligo- and polycarbonic acids and other light-cure dental resins. The ionomer glass used contains fluoride and zinc ions. Filler content: 72% by weight = 55% by vol. (0.02–33 µm).

Packaging

Packaging	
3 Saferinge containing 0.33 g paste	REF 220091
20 Saferinge containing 0.33 g paste	REF 209538
2 syringes containing 1.5 g paste, 20 Luer-Lock-Tips	REF 213002
6 syringes containing 1.5 g paste, 40 Luer-Lock-Tips	REF 213001

Accessories	
60 Luer-Lock-Tips	REF 211477

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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