# LuxaBite

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English

### Product description

LuxaBite is a bis-acrylic based self-mixing bite registration material. It offers barely discernible resistance to the jaw occlusion of the patient. As such, a reliable bite registration of the occlusion is guaranteed. A detailed recording of the occlusal relief and the high end hardness simplify the precise positioning of working and counter-bite models.

### Indications

Bite registration

### Basic safety instructions

- ▶ For dental use only!
- ▶ Keep out of reach of children!
- ▶ Avoid contact with the skin! In the event of inadvertent skin contact, wash the affected area thoroughly with soap and water.
- ▶ Avoid contact with the eyes! In the event of inadvertent contact with the eyes, rinse thoroughly with plenty of water immediately and consult a physician if necessary.
- ▶ Do not use LuxaBite in the case of pre-existing hypersensitivity to acrylates.
- ▶ Do not apply LuxaBite directly to the oral mucosa in tooth-free areas since heat may develop.

### Interactions

- ▶ Avoid contact with retraction solutions, because in some instances residues of retraction solutions can increase the setting time. Rinse residues of retraction solutions thoroughly with water prior to the bite registration.

### Using the Automix cartridge

#### Insert cartridge

1. Push up the lever [Fig. 1, Pos. A] at the back of the Automix-Dispenser and pull back the slide [Fig. 1, Pos. B] completely.
2. Lift up the plastic lever [Fig. 1, Pos. C], insert the cartridge and lock it using the plastic lever [Fig. 1, Pos. C].

**Note:** Make sure that the notches on the cartridge line up with those on the Automix-Dispenser.

#### Attach mixing tip

1. Remove the cartridge cap [Fig. 1, Pos. D] or used mixing tip by turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [Fig. 1, Pos. E].

**Note:** In order to attain an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are also available as refill packs.

**Note:** Make sure that the notches on the mixing tip line up with those on the cartridge.

3. Turn the mixing tip clockwise 90° to lock it.

#### Applying the material

- ▶ When using a cartridge for the first time, squeeze out a pea-sized amount of material and discard it.
- ▶ The material is mixed in the mixing tip by pulling the trigger [Fig. 1, Pos. F] of the Automix-Dispenser and it is then applied directly.

**Note:** After the application, leave the used mixing tip on the cartridge to cap unused material!

#### Remove cartridge

1. Push up the lever [Fig. 1, Pos. A] at the back of the Automix-Dispenser and pull back the slide [Fig. 1, Pos. B] completely.
2. Lift up the plastic lever [Fig. 1, Pos. C] and remove the cartridge.

### Instructions for use

- ▶ For the bite registration of complete dental arches, use LuxaBite in segments (3 – 4 cm).
- ▶ For the individual times refer to the “Technical Data” section.

### Recommended use

#### Habitual bite registration

1. Lightly isolate any existing acrylic structures with Vaseline or a similar separating agent. Block out any existing bridge pontics using a suitable substance, e.g. wax.
2. Apply LuxaBite to the dry dental row. When doing so the length of the material may be 3 to 4 cm. Application should take place to prepared areas first.
3. Allow patient to adopt the desired occlusion.
4. In order to ensure the complete bite registration of the occlusal surfaces, apply additional material under the buccal cusps of the maxillary teeth.
5. Once the material has cured for 45 to 60 s, have the patient open their bite slightly and then close it once more. This will ease the subsequent removal of the register.
6. Once the material has cured for a further 90 s, remove the material with gloved fingers or using a suitable instrument.
7. In order to avoid positioning errors, reduce the register to the occlusal contact surfaces. To do so use a suitable instrument, e.g. an acrylic bur.
8. Articulate models with the aid of the bite register.

#### Central bite registration using anterior tooth jig

1. Use a suitable material, e.g. composite or silicone, for the vertical opening. When doing so, position the material at the anterior mandible and have the patient assume the desired occlusion. The dental rows are not fully closed, in order that space remains for the application of the bite registration material.
2. The procedure now continues as described under points 2 to 8 for habitual bite registration.

### Disinfection of the bite registration

The bite registration can be disinfected using standard disinfectant solutions. The use of products containing alcohol is recommended because these additionally remove the oxygen inhibition layer.

### Storing the bite registration

- ▶ Protect against exposure to direct sunlight.
- ▶ Store at room temperature (15 – 25 °C/59 – 77 °F) in a dry place.

The bite registration can be stored without time limitation.

### Technical Data

Total working time	0:45 min from starting mixing
Elastic phase	2:00 to 2:30 min from starting mixing
Final curing	4:00 min after removal
Barcol hardness	25 after 1 h

### Disinfection of the Automix-Dispenser

The Automix-Dispenser is reusable after disinfection, although it should be replaced if visible signs of wear appear.

Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfecting solutions listed by the RKI (Robert Koch Institute). The Automix-Dispenser is also autoclavable.

### Storage and shelf life

- ▶ Store at room temperature (15 – 25 °C/59 – 77 °F) in a dry place.
- ▶ Do not store in the fridge or freezer.
- ▶ Protect against exposure to direct sunlight.
- ▶ After the application, leave the used tip on the cartridge to cap unused material.
- ▶ Do not use after the expiration date.

## Ingredients

Glass filler material in a matrix of multifunctional methacrylates. Free from methyl methacrylate and peroxides.

## Available product sizes

1 cartridge @ 50 ml paste, 15 Automix-Tips	REF 110560
5 cartridges @ 50 ml paste, 75 Automix-Tips	REF 110558

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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