LuxaPrint

Cast

DMG Chemisch-Pharmazeutische Fabrik GmbH
 Elbgaustraße 248 · 22547 Hamburg · Germany
 www.dmg-dental.com
 092040/#1_2016-12

DE	Gebrauchsinformation	IT	Istruzioni per l'uso
EN	Instructions for use		

Gebrauchsinformation

Deutsch

Produktbeschreibung

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English

Product description

(Meth) acrylate-based light-curing resin exclusively for use in the dental laboratory. For DLP / STL printers that work at a wavelength of 405 nm or 385 nm.

Purpose

Manufacture of molds that are combustible without any residue (e.g. model casting frameworks, crowns and bridges) using 3D printing.

Restrictions on use

- Intended exclusively for the creation of form works and structures that are combustible without any residue in the dental sector.
- Not a medical device according to EG “Medical Devices” Directive 93/42/EEC.
- May not be brought into direct contact with patients.

Basic safety information

- Improper use and failure to observe information may lead to impaired quality.
- Observe the safety data sheet!
- Keep out of the reach of children!

- Wear safety gloves (nitrile gloves), protective clothing and safety goggles while processing.
- Avoid eye contact! In the event of the liquid material accidentally coming into contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Avoid direct skin contact with the liquid material and components prior to post-curing.

Side effects

- Irritates the eyes and skin (sensitization possible).
- Do not use the polymerized material if you are allergic to any of the ingredients or, as may occur in rare cases, in the event of contact allergies.
- Avoid skin contact with non-polymerized material and the inhalation of monomer vapors. In rare cases, allergic reactions to LuxaPrint Cast components may occur.

Notes for use

For optimum properties on the end product:

- Ensure that the light-curing unit is in a proper condition.
- Leave objects to cure completely.
- Maintain a processing temperature of $23^{\circ}\text{C} \pm 2^{\circ}\text{C}$ ($73^{\circ}\text{F} \pm 2^{\circ}\text{F}$).
- Ensure clean working. Residues of dirt on the machine can cause defects on the printed object.
- Observe the machine and material parameters.

Recommended use

Note: Shake the bottle for at least 1 min before use.

Reworking the printed object

Work steps	Time [min]
Recommended drip-off time after raising the platform	10 min
Note: Reworking should take place immediately after the construction process when possible.	
1. Loosen the printed object from the construction panel. Carefully loosen the support structures and support.	
Optional: Pre-clean the printed object using ethanol (approx. 96%) and a brush.	3 min
Optional: Clean openings, holes and gap areas with compressed air.	
2. Clean the printed object in a separate container using fresh ethanol in an ultrasonic bath.	3 min
Optional: Heat the printed object up to approx. 40°C (104°F) in the oven to remove solvent residue from the cleaning process	30 min

Note: Burn the printed object on a model to avoid separation of the support structures.

Burning the printed object

Burning takes place (e.g. in the Otoflash G171 flash-curing unit) in a protective gas atmosphere (nitrogen). Observe the manufacturer's instructions for the unit.

Examples of light-curing units / times

Light-curing unit	Light-curing time*	Total output [W]
Otoflash G171 Xenon flash-curing unit	2 x 2,000 flashes Turn the printed object in between!	Two 100 W flash tubes arranged at the bottom
Blue UV-A	10 min	71
Heraflash	2 x 180 s Turn the printed object in between!	

* The information on the light-curing times is based on the manufacture of lost molds from DMG material with a DLP 405 nm 3D printer.

Curing temperatures

Please observe the manufacturer's instructions for the embedding compound used. A maintenance level of 850°C must be maintained for roughly 30–45 min (depending on the size of the embedded objects), regardless of the final casting temperature, so that the resin cures completely.

Composition

Mixture of acrylate and methacrylate resins, photoinitiators, colorants and additives.

Classification

LuxaPrint Cast is not a medical device according to the EG "Medical Devices" Directive 93/42/EEC. The product may not be brought into direct contact with patients.

Storage and shelf life

- Dry at room temperature ($15^{\circ}\text{--}25^{\circ}\text{C}$ / $59^{\circ}\text{--}77^{\circ}\text{F}$) and store in a location that is protected against light!
- Even low exposure to light can trigger polymerization:

Note: While handling (e.g. transferring), only expose LuxaPrint Cast to daylight / ambient light for as short a time as possible, since exposure to light for too long a time can lead to undesirable curing of the material.

- Always keep the container tightly sealed; carefully seal immediately after every instance of use.
- Do not use after the expiration date!

Packaging

1 kg liquid	REF 170203
0.5 kg liquid	REF 170202

Warranty

DMG is certified according to DIN EN ISO 13485 and operates an elaborate quality assurance system. It is therefore guaranteed that this product is free from material and manufacturing defects and is thus of a perfect quality.

Our user recommendations are based on careful determination of recommended values from our test laboratory. To achieve the guaranteed level of product quality and usability, the specified process steps must be strictly complied with. Once the product has been processed, DMG assumes no liability and provides no guarantee with respect to the saleability or suitability of the product for a particular purpose. The user is responsible for the product's application and intended use.

If defects occur on the product within the warranty period, the user is only entitled to a product replacement. DMG assumes no liability for losses or damages caused by this product – irrespective of whether they are direct, indirect, special, collateral or consequential damages – whatever the legal basis. DMG does, however, assume liability for direct material damages to the product that are attributable to intent or gross negligence on the part of its legal representatives or executive employees and for personal injuries in accordance with statutory provisions. Any liability for the product and consequential damages arising from use of the same is excluded if the user has failed to observe the specified process steps.

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