LuxaPrint

Ortho

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Product description

(Meth) acrylate-based light-curing resin exclusively for use in the dental laboratory. For DLP / STL printers that work at a wavelength of 405 nm or 385 nm.

Indications

Generative manufacture of drilling templates and orthodontic equipment using 3D printing.

Contraindications

Only insert LuxaPrint Ortho intra-orally in a fully polymerized state.

Basic safety information

- Improper use and failure to observe information may lead to impaired quality.
- Observe the safety data sheet!
- Keep out of the reach of children!
- Wear safety gloves (nitrile gloves), protective clothing and safety goggles while processing.
- Avoid eye contact! In the event of the liquid material accidentally coming into contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Avoid direct skin contact with the liquid material and components prior to post-curing.
- Biocompatibility is only guaranteed in the event of proper use (polymerization)!

Side effects

- Irritates the eyes and skin (sensitization possible).
- Do not use the polymerized material if you are allergic to any of the ingredients or, as may occur in rare cases, in the event of contact allergies.
- Avoid skin contact with non-polymerized material and the inhalation of monomer vapors. In rare cases, allergic reactions to LuxaPrint Ortho components may occur.

Notes for use

For optimum properties on the end product:

- Ensure that the light-curing unit is in a proper condition.
- Leave objects to cure completely.
- Maintain a processing temperature of $23^{\circ}\text{C} \pm 2^{\circ}\text{C}$ ($73^{\circ}\text{F} \pm 2^{\circ}\text{F}$).
- Ensure clean working. Residues of dirt on the machine can cause defects on the printed object.
- Observe the machine and material parameters.
- Do not use heat-based disinfection or sterilization methods. Otherwise, the possibility of the workpiece deforming cannot be eliminated.
- Minimum material thickness for structures: 1.5 mm
- Maximum curing depth during burning:
For solid objects and light-curing on both sides, the material can be up to 7 mm thick (at a curing depth of 3.5 mm).
- Select material variants to suit the wavelength of the printer used (385 nm or 405 nm).

Recommended use

Note: Shake the bottle for at least 1 min before use.

Reworking the printed object

Work steps	Time [min]
Recommended drip-off time after raising the platform	10 min
Note: Reworking should take place immediately after the construction process when possible.	
1. Loosen the printed object from the construction panel. Carefully loosen the support structures and support.	
Optional: Pre-clean the printed object using ethanol (approx. 96%) and a brush.	3 min
Optional: Clean openings, holes and gap areas with compressed air.	
2. Clean the printed object in a separate container using fresh ethanol in an ultrasonic bath.	3 min
Optional: Heat the printed object up to approx. 40°C (104°F) in the oven to remove solvent residue from the cleaning process	30 min
3. Burning: Burn the printed object under a protective gas atmosphere (nitrogen) and without heating (e.g. in the Otoflash G171 Xenon flash-curing unit) to prevent yellowing of the components. Observe the manufacturer's instructions for the unit. Note: Burn the printed object on a model to avoid separation of the support structures.	
4. Processing the surface (e.g. painting)	

Examples of light-curing units / times

Light-curing unit	Light-curing time*	Total output [W]
Otoflash G171 Xenon flash-curing unit	2 x 2,000 flashes Turn the printed object in between!	Two 100 W flash tubes arranged at the bottom
Blue UV-A	5 min	71

Heraflash	2 x 180 s Turn the printed object in between!	
-----------	--	--

Composition

Mixture of acrylate and methacrylate resins, photoinitiators and additives.

Classification

Class 1

Storage and shelf life

- Dry at room temperature (15°–25° C / 59°–77° F) and store in a location that is protected against light!
- Even low exposure to light can trigger polymerization:

Note: While handling (e.g. transferring), only expose LuxaPrint Ortho to daylight / ambient light for as short a time as possible, since exposure to light for too long a time can lead to undesirable curing of the material.

- Always keep the container tightly sealed; carefully seal immediately after every instance of use.
- Do not use after the expiration date!

Packaging

385 nm	
1 kg liquid	REF 170207
0.5 kg liquid	REF 170206

405 nm	
1 kg liquid	REF 170209
0.5 kg liquid	REF 170208

Warranty

DMG is certified according to DIN EN ISO 13485 and operates an elaborate quality assurance system. It is therefore guaranteed that this product is free from material and manufacturing defects and is thus of a perfect quality.

Our user recommendations are based on careful determination of recommended values from our test laboratory. To achieve the guaranteed level of product quality and usability, the specified process steps must be strictly complied with. Once the product has been processed, DMG assumes no liability and provides no guarantee with respect to the saleability or suitability of the product for a particular purpose. The user is responsible for the product's application and intended use.

If defects occur on the product within the warranty period, the user is only entitled to a product replacement. DMG assumes no liability for losses or damages caused by this product – irrespective of whether they are direct, indirect, special, collateral or consequential damages – whatever the legal basis. DMG does, however, assume liability for direct material damages to the product that are attributable to intent or gross negligence on the part of its legal representatives or executive employees and for personal injuries in accordance with statutory provisions. Any liability for the product and consequential damages arising from use of the same is excluded if the user has failed to observe the specified process steps.

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